

CONFIDENTIAL B

The Mediatek logo consists of the word "MEDIATEK" in white, uppercase, sans-serif font, centered within an orange parallelogram shape that is wider on the left and tapers to the right.

**MEDIATEK**

# MT2511 Health Module Programming Guide

**2017.04.19**

# Outline

- Overview
- Introduction
  - Main Class
  - Send Command Flow
  - Receive Data Flow
  - Data Format
- Appendix

# Overview

# Overview

- **MT2511 Health Module**


- bases on *MT2511Controller*@Smart Device, could receive feature data(Heart Rate/Blood Pressure/PWTT) and signal raw data (EKG/PPG1/PPG1-512hz) from device.

(Note: receive signal raw data only in APK SPP mode.)

- It provides the following features:

- Add record – Transfer personal profiles into device
- HR Measure – measure heart rate
- BP measure
  - General Mode – Measure personal blood pressure
  - Personal Mode – Measure bp with saved personal model, which created from calibration process.
- Multiple connect – Add reference device for heart rate comparison

# Architecture



# Introduction

# Package Structure

- *Package Structure*
  - DemoKit project include:
    - `iot`: Data Transfer & Bluetooth features  
BTData & DataParser: Data format module
    - `mwcdemo`: MT2511 UI features

# Main Class (1/3)

- *MT2511Controller*
  - Inherit from Wearable SDK *Controller*
  - Receiver: *health\_receiver*.
  - Implement *onConnectionStateChange*, *onReceive* to receive and post MT2511 *ReceiveData* by using *RxBus*.
  - Add/Remove *MT2511Controller* in *MainService onCreate/onDestroy*.
  - *Handshake* with device to Determine device whether or not support MT2511 health feature.
  - Observe *RequestWriteToDevice* by using *RxBus* and use *Controller send* API to *sendCommand*.



# Main Class (2/3)

## ■ *AdapterDevice*

- Inherit from DemoKit/iot *Device*.
- Observe *ReceiveData* from *MT2511Controller* and using *BTDataParser* to parse it, then post *BTBaseData* (HR/BP/EKG etc.) to UI module by using *RxBus*.
- Provide *writeToDevice* API for UI module, *writeToDevice* could post *RequestWriteToDevice* to *MT2511Controller* and let *MT2511Controller* to send command.
- *Connect/Disconnect* only post event by using *RxBus*.
- Using *AdapterDeviceFactory* to crease *AdapterDevice*.

```
.sBTDevice = new AdapterDevice(MContext.getInstance().getApplication(),  
    | ..... new BTDataParser());
```

# Main Class (3/3)

## ■ *RxBus*

- Data Bus implement with *RxJava*, could *post* java Object <T> to observer (call *toObservable* <T> class).
- For more *RxJava* info, please refer to Appendix *RxJava*.

- Sample Code:

Subject Poster — *DataParser/BTDataParser postData*

Observer — *MeasureFragment initView*


```
protected void postData(byte[] data) {  
    try {  
        BaseData baseData = parseData(data);  
        if (isRememberData(baseData)) {  
            RxBus.getInstance().postAndRemember(baseData);  
        } else {  
            RxBus.getInstance().post(baseData);  
        }  
    }  
}
```

```
_subscriptions.add(RxBus.getInstance()  
    .toObservable(HeartRateData.class)  
    .observeOn(AndroidSchedulers.mainThread())  
    .subscribe(new Action1<HeartRateData>() {  
        @Override public void call(HeartRateData heartRateData) {  
            int hr = heartRateData.get(HeartRateData.FIELD_BPM);  
            String bmp_view_value = String.valueOf(hr);  
            mTxtHeartRate1.setText(bmp_view_value);  
        }  
    }  
));
```

# Send Command Flow (1/3)

- Send Command Flow – Send Person Profile

Add Record  
(send person profile  
info to device)




The screenshot shows a mobile application interface with a status bar at the top displaying the time 13:49, network speed 0.42K/s, and various icons. The app's title bar is brown and contains a back arrow, the text "BP Measure", and a plus icon. The main content area is titled "Add Record" in orange. It contains several input fields: "User ID" with the value "Jack", "Age" with the value "29", "Height (cm)" with the value "173", and "Weight (kg)" with the value "68". There are radio buttons for "Male" (selected) and "Female". At the bottom of the form are "CANCEL" and "OK" buttons. Below the form is a brown "START" button.

# Send Command Flow (2/3)

- Send Command Flow – Send Person Profile


*sendCommand*  
(Base Class *Device*  
method)

Run in  
*CommandOnSubscribe*




# Send Command Flow (3/3)

- Send Command Flow – Send Person Profile




# Receive Data Flow (1/3)

- Receive Data Flow – Receive HR bmp




# Receive Data Flow (2/3)

- Receive Data Flow – Receive HR bmp



# Receive Data Flow (3/3)

- Receive Data Flow – Receive HR bmp





# Data Format (1/8)

- **Command** (APK -> Device)
  - All command type extend from [BaseCommand].
  - Class Hierarchy




- Sample Code:
  - *DemoKit\src\main\java\com\mediatek\iot\command*
  - *CustomUtils*

# Data Format (2/8)

## ■ Data (Device -> APK)

- All data type extend from [BaseData].
- Data Size 64 bytes = 16 \* int (4 bytes)
- Class Hierarchy



## – Sample Code:

- `DemoKit\src\main\java\com\mediatek\iot\data`
- `BTDataParser`
- User can use `BTBaseData.get(int index)` method, get the filed as user need.

```
int hr = heartRateData.get(HeartRateData.FIELD_BPM);
```

# Data Format (3/8)

## ■ EKGDData

- A int array with length = 12 of EKG raw data
- The EKG data are recorded in 512 hz.

EKG
[0] magic (54321)
[1] sensor type ( sensor type is 5 )
[2] sequence
[3-14] 12 ekg data
[15] reserved (12345)

```
public int[] getEkgData() {  
    ... int[] data = new int[12];  
    ... System.arraycopy(rawData, 3, data, 0, 12);  
    ... return data;  
}
```

## – Sample Log:

54321,5,111,7403228,7401915,7405257,7408628,7407863,7408847,7410483,7413195,7413959,7415612,7415998,7413969,12345

54321,5,112,7415451,7418458,7418712,7419345,7419505,7417862,7417064,7419372,7419947,7420994,7422060,7420012,12345

# Data Format (4/8)

## ■ PPG1Data

- A int array with length = 12 of PPG1 raw data
- The PPG1 data are recorded in 125 hz.

PPG1
[0] magic (54321)
[1] sensor type ( sensor type is 9 )
[2] sequence
[3-14] 6 ppg + 6 amb
[15] reserved (12345)

### – Sample Log:

54321,9,766,8356541,8363026,8356541,8363004,8356541,8363016,8356541,8363018,8357172,8363010,8359011,8363056,12345,  
54321,9,767,8360529,8363029,8361763,8363007,8362794,8363034,8363612,8363030,8364296,8363062,8364873,8363123,12345,  
54321,9,768,8365327,8363074,8365703,8363078,8365977,8363055,8366206,8363049,8366450,8363071,8366610,8363037,12345,

# Data Format (5/8)

## ■ PPG1\_512Data

- A int array with length = 12 of PPG1 raw data
- The PPG1 data are recorded in 512 hz.

PPG1_512HZ
[0] magic (54321)
[1] sensor type ( sensor type is 12 )
[2] sequence
[3-14] 12 PPG Data
[15] reserved (12345)

- Sample Log:

54321,12,1331,6018,5995,6005,6038,6052,6033,5997,6004,6044,6052,6032,5985,12345,  
54321,12,1332,5991,6027,6046,6037,5996,6001,6029,6056,6041,6007,6004,6037,12345

# Data Format (6/8)

## ■ HeartRateData

HR
[0] magic (54321)
[1] sensor type ( sensor type is 22 )
[2] sequence
[3] bpm
[4] status
[5] timestamp
[6-15] reserved (12345)

```
public static final int FIELD_BPM . . . . . = 3;  
public static final int FIELD_STATUS . . . = 4;  
public static final int FIELD_TIMESTAMP = 5;
```

// HR Confidence Level = status & 0x000000FF

— Sample:

— *MeasureFragment toObservable(HeartRateData.class)*

54321,22,105,92,1228341247,-43075096,12345,12345,12345,12345,12345,12345  
,12345,12345,12345,12345

# Data Format (7/8)

## ■ BloodPressureData

Blood Pressure
[0] magic (54321)
[1] sensor type ( sensor type is 24 )
[2] sequence
[3] sbp
[4] dbp
[5] HR bpm
[6] status
[7] timestamp
[8-15] reserved (12345)

```
public static final int FIELD_SBP.....= 3;  
public static final int FIELD_DBP.....= 4;  
public static final int FIELD_HR_BPM....= 5;
```

– Sample:

– *BPMeasureFragment showBPMeasureResult*

54321,24,1,109,76,62,0,0,12345,12345,12345,12345,12345,12345,12345,12345

# Data Format (8/8)

## ■ PWTTData

- FEATURE\_TYPE = 0:  
int array represent the **PWTT interval**  
of the blood pressure measure result
- FEATURE\_TYPE = 3:  
int array represent the **personal model**  
of the blood pressure calibration process

- Sample:

- *BPMeasureFragment.toObservable(PWTTData.class)*

54321,8001,1,0,9,0,192,162,180,204,56,153,210,53,259,12345

54321,8001,2,0,7,0,254,59,150,153,215,184,208,53,259,12345

54321,8001,3,3,6,0,109,76,62,182,324,0,208,53,259,12345

PWTT (pulse wave transit time)	
[0]	magic (54321)
[1]	sensor type (sensor type is 8001 )
[2]	sequence
[3]	feature type
[4]	status
[5]	timestamp
[6..N]	feature value
(N<15)	
[N+1..15]	reserved (12345)



# Appendix

# Dependency

## ■ RxJava

- A library for composing asynchronous and event-based programs using observable sequences for the Java VM.
- Web: [Official Web](#) in Github [How-To-Use-RxJava](#)
- Build.gradle :

```
// RxJava - .composing asynchronous and event-based programs using observable sequences  
compile 'io.reactivex:rxjava:1.1.0'  
compile 'io.reactivex:rxandroid:1.1.0'  
compile 'com.jakewharton.rxbinding:rxbinding:0.3.0'
```

- Sample Code
  - *RxBus/DataParser/AdaperDevice* Java File
  - *MeasureFragment initView()*

# SNR Computation

- **SNR** is comparison result of the level of signal to the level to noise.
- SNR Class:  
DemoKit\demoKit\com\mediatek\mwcdemo\snr\SNRResult.java
- Sample Code:  
DemoKit\demoKit\com\mediatek\mwcdemo\fragments\BPMeasureFragment.java

```
// insert high pass filtered input to the lib when raw data received
```

```
double s_hpf = mECGFilterService.filter(data_ecg_512_mv);
```

```
mECGSNR.inputHPFSignal(s_hpf);
```

```
...
```

```
// when the measurement is finished, compute the low pass signal array with convolution function.
```

```
//compute the SNR result afterwards.
```

```
hpfSignalList = mECGSNR.getHPFSignal();
```

```
allLPFSingalList = mECGFilterService.conv(hpfSignalList);
```

```
mECGSNR.computeECG512SNR40(allLPFSingalList);
```



*everyday genius*